

MCC Calculator Workflow Procedures: REG / CBF / ASL

Contents

- REG Tool – Elastic Registration
- CBF Measurement
- ASL Measurement

Installation & Setup Prerequisites


- MATLAB installed and properly licensed
- Access to sample datasets / image sets used by MCC tools
- Sufficient disk space for importing and exporting image sets

MATLAB Requirements


- MATLAB R2020b or later
- Image Processing Toolbox
- Signal Processing Toolbox (for CBF analysis)
- Parallel Computing Toolbox (recommended for REG tool)

Procedure - MCC_Calculator_REG - Elastic Registration

1. Launch **RegCalculator.mlapp**.
2. Click **Import** to load the image set.




3. Browse and select the **target sample data folder**.
4. The **reference image** and **threshold** will be calculated automatically. Click **Register** to start elastic registration.




5. Wait for the registration to complete (typically **5 – 8 minutes**, depending on dataset size and PC performance).

6. Review the registered results in the **left preview panel**.




7. Click **Export** to save the registered image set for further processing.




Procedures - MCC_Calculator_4 - CBF Measurement


1. Launch **MccCalculator.mlapp**.
2. Click **Load** to import the image set.




3. Browse to and select the target sample data folder.
4. Go to the **CBF** tab (top-right corner of the interface).



5. Using dynamic observation of the image sequence, manually define:
- (1) the **cilia region (ROI)**, and
 - (2) a **background region (ROI)** for reference.




6. Click **Calculate** to run the CBF computation.
7. The calculated CBF result will be displayed in the textbox at the bottom-right of the interface.




Procedures - MCC_Calculator_4 - ASL Measurement

1. Launch **ASLCalculator.mlapp**.
2. Click **Import** to load the image set.



3. Browse to and select the target sample data folder. Choose the node to load original image set from **left panel**.
4. Click **Segment** on the ribbon bar, then manually segment (draw) the two layers in the image:
 - (1) the **upper surface layer**, and
 - (2) the **epithelium layer**.



5. Click **Process** to start the ASL measurement.
6. Wait about **2 minutes** for layer segmentation/processing to complete.
7. Review the calculated reference ASL under the “**ASL**” node in the left panel. The ASL thickness result will be displayed in the **textbox at the bottom-right** of the interface.

